

TECHNICAL DATA SHEET

Safetouch® Vitals™

Wash glove – Needle punch
70g/m²

Kolmi®



FEATURES AND BENEFITS

- Absorbent, soft and very resistant washing glove
- Ultrasonic sealing
- Compact & hygienic vacuum polybag

RECOMMENDED USE

- Used for body hygiene

PRODUCT DESCRIPTION

Product name: Washglove – Needle punch
Product type: Single use, non-sterile
Material: Needle punch nonwoven 100% Polyester
Ultrasonic sealing
Colour: White
Packaging: 70 or 50 units in vacuum polybag
Country of origin: China



Article	Color	Dimensions (cm)	Unit weight (g)	Packaging
G01421	White	14 x 21	4,1	20 bags x 70 units
G01521		15 x 21	4,3	20 bags x 50 units
G01623		16 x 23	5,5	

PRODUCT PERFORMANCES

Test	Standard		Average	Unit	Tolérance
Grammage	ISO 9073-1		70	g/m ²	± 10%
Thickness	ISO 9073-2		0,33	mm	± 10%
Tensile Strength	ISO 9073-3	Wet MD:	67	N/5cm	± 10%
		Wet CD:	44	N/5cm	± 10%
Elongation	ISO 9073-3	Wet MD:	56	%	± 10%
		Wet CD:	70	%	± 10%
Absorption	ISO 9073-6		500	%	± 20%

STORAGE CONDITIONS

Normal conditions of conservation & storage: must be stored in a cool or room temperature environment and in a dry atmosphere.

PACKAGING MARKING



* Do not use if packaging is damaged

LOGISTICAL INFORMATION

OUTER CASE SPECIFICATIONS

ITEM	SIZE (mm)	WEIGHT (kg)	QTE/ PALLET	OUTER CASE EAN	QTE / OUTER CASE	POLYBAG EAN
G01421	300 x 220 x 430	6	36 boxes 6 layers of 6	3 662 036 021 243	20 Polybags of 70 units	3 662 036 021 236
G01521	320 x 220 x 320	4,7	60 boxes 5 layers of 12	3 662 036 021 229	20 Polybags of 50 units	3 662 036 021 212
G01623	340 x 240 x 320	5,5	60 boxes 10 layers of 6	3 662 036 019 738		3 662 036 019 745