

# Microlab 615 Basic Dual Syringe Diluter

Specification Sheet (Part/REF # ML615-DIL)



The Microlab 600's Dual Syringe Diluter configuration uses two syringes to create up to a 1:50,000 dilutions in a single step, drastically reducing preparation time and wasted reagent. The diluent washes the tubing between each sample, minimizing carryover.

## Product Specifications

<b>Power Rating</b>	24 VDC, 2.5A
<b>Power Requirements</b>	100-240V 1.5A max 50-60 Hz
<b>Program Memory</b>	One program stored in memory while power is on
<b>Syringe Drive Mechanism</b>	Stepper motor driven high precision lead screw with encoder
<b>Controller</b>	Basic
<b>Number of Syringe Drives</b>	Two
<b>Probe</b>	Concorde CT Hand Probe

Specifications are subject to change without notice

**Spec. Version B**

<https://www.hamiltoncompany.com/laboratory-products/microlab-600/ml615-dil>