

TECHNICAL DATA SHEET

Part Number	49199-50ML-HON
Description	Acetic acid

Technical specifications :

- ✓ acetic acid, eluent additive for lc-ms 50ml
- ✓ CAS N°64-19-7
- ✓ Molecular weight 60.05 g/mol
- ✓ Molecular formula C₂H₄O₂
- ✓ Density 1.05 g/cm³ (20 °C)
- ✓ Flash point 40 °C
- ✓ Material form liquid

- ✓ Storage condition Room Temperature
- ✓ Package type Glass Bottle
- ✓ Package size 50mL

- ✓ GHS pictograms GHS02 ; GHS05
- ✓ H Statements H314 ; H226
- ✓ P Statements P260 ; P280 ; P284

Sales unit :

50mL Acetic acid
