



PolyPress™ Tissue Embedding Tampers



← Long handle for easy gripping



Available in two sizes

PolyPress™ Tissue Embedding Tampers aid in the orientation of samples during embedding. There are three options: small, large and a combination pack.

Features:

- Easy to clean and reduces build up
- Metal keeps paraffin from sticking to fingers
- 2 Convenient sizes: large - 15mm x 15mm & small - 10mm x 10mm
- No need for forceps

Catalog #25960

Size: small

Catalog #25964

Size: large

Catalog #26013

Size: small/large combination