

## TECHNICAL DATASHEET

Ref. DD **498401** - code fabricant : 3500922

Mop holder Clino® Synchro PRO CR 40 cm

Mop frame for Cleanroom

Individualization on request!

**NEW**



### Product description

Newly developed, autoclavable mop holder with robust locking mechanism made of high-performance plastic for the cleanroom.

- Hygienic design
- Extremely robust closing mechanism
- Pre-tensioning of the wings for optimal pressure distribution
- Easy mop pick up
- Very flat guide possible

**MADE IN GERMANY**

<b>ISO</b>	3	4	5	6	7	8	9
<b>GMP</b>			A/B		C	D	

### Product specifications

Article name   Article number	Mop holder Clino® Synchro PRO CR 40 cm   3500922
Bemaßung   Weight	40 x 10 cm   360 g
Opening symmetry	Symmetric
Wing material	Polyamide 6
Opening mechanism material	POM
Connection mechanism	Clino® Connect Lock
Connection Mechanism material	Polyamid 6 (Heat-stable special plastic)
Sterility	Autoclavable when open
Packing unit   Selling unit	1 Pieces   1 PU
Quantity in carton	20 Pieces
	7200 g
Pallet	120 x 80 x 220 cm
Packaging note	The size of the package, as well as the quantity in the box may vary depending on the order quantity and the volume of other ordered products.
Customs tariff number	96039010
Country of origin	Germany
EAN-Number	4065237110584
Compatible accessories	Our sales department will be happy to advise you on suitable accessories

### Compatible methods

- MopScoop
- MopFloat
- EasyMop
- Flat wringer Ringo

### Cleaning area



### Industry Sectors



### Hygiene



Diese Klassifizierung ist eine unverbindliche Empfehlung. Der Anwender muss die individuelle Eignung und Verträglichkeit für seinen individuellen Anwendungsbereich prüfen. Der Zustand des Materials ist regelmäßig zu überprüfen, beschädigte Produkte sind umgehend auszutauschen. Alle Angaben ohne Gewähr. Änderungen und Abweichungen vorbehalten.