


Ref. DD : 498402 - Code fabricant : 2009461

Mop holder stainless steel inverse with Clino® Connect Lock

Mop frame for Cleanroom



Product description

Autoclavable mop holder made of sturdy and durable stainless steel and high performance plastic specially designed for disposable mops.

- o Autoclavable
- o Stainless steel plate at the bottom
- o Not suitable for flat wringer

MADE IN EUROPE

ISO	3	4	5	6	7	8	9
GMP	A/B	C	D				

Product specifications

Article name | Article number
Dimensions | Weight
Opening symmetry
Wing material
Opening mechanism material
Connection mechanism
Connection Mechanism material
Sterility
Packing unit | Selling unit
Quantity in carton
Carton weight
Pallet
Packaging note

Mop holder stainless steel inverse with Clino® Connect Lock | 2009461
40 x 10 cm | 690 g
Asymmetric
Stainless steel 1.4401
Stainless steel 1.4301
Clino® Connect Lock
Polyamid 6 (Heat-stable special plastic)
Autoclavable
1 Pieces | 1 PU
20 Pieces
13800 g
120 x 80 x 220 cm

The size of the package, as well as the quantity in the box may vary depending on the order quantity and the volume of other ordered products.

Customs tariff number
Country of origin
EAN-Number
Compatible accessories

96039010
England
4065237003688
Our sales department will be happy to advise you on suitable accessories

Compatible methods Cleaning area

- o MopScoop
- o MopFloat
- o EasyMop
- o Flachpresse Ringo



Industry Sectors

Hygiene

Diese Klassifizierung ist eine unverbindliche Empfehlung. Der Anwender muss die individuelle Eignung und Verträglichkeit für seinen individuellen Anwendungsbereich prüfen. Der Zustand des Materials ist regelmäßig zu überprüfen, beschädigte Produkte sind umgehend auszutauschen. Alle Angaben ohne Gewähr. Änderungen und Abweichungen vorbehalten.

PFENNIG Reinigungstechnik GmbH | Heubachstraße 1 | 87471 Durach | Germany
www.pfennig-reinigungstechnik.com | info@pfennig-reinigungstechnik.com | +49 (0) 831 561 22 0

Status: 25.10.2023

