



AmpliPur Expert Tips

Part Number	F174401
-------------	---------

Specification

Description	AmpliPur Expert Tips, 100-1000 µL, graduated, filter tips
Dimensions	Tip length: 88 mm Inner diameter tip point: 0.9 mm Outside tip collar diameter: 8.7 mm Rack box dimensions: 122.5 x 90.5 x 104 mm
Volume Range	100-1000 µL
Nominal Volume of Compatible Pipettes	1000 µL
Material (Plastic)	<ul style="list-style-type: none">• Tips: Polypropylene (PP)• Filter: Polyethylene (PE)• Rack: Polypropylene (PP)• Heat shrink: Plastic Optical Fibre (POF)
Disposable	<ul style="list-style-type: none">• Tips: yes• Rack: reusable
Tip Color	Translucent
Sterilization	N/A
Quality Control	Material control Product inspection: control of characteristics according Gilson's specifications

Characteristics

Characteristic	The optimized geometry and material (virgin polypropylene) ensures full recovery of the aspirated liquid for consistent results
Temperature	Maximum temperature 140°C
Autoclavability	Not applicable on filter tip
Expiration Date/Shelf Life	3 years
Traceability/Quality	<ul style="list-style-type: none">• Complete traceability from production to the user• Certificate of conformity• Quality statement with protocols and detection thresholds for each batch
Other Information	<ul style="list-style-type: none">• Compliant to ISO 8655-2 maximum permissible errors• Free of RNase and DNase• Compliant to the use of hazardous substances directive 2015/863/EU• Graduation mark at 50 µL, 100 µL, 200 µL, and 500 µL

Packaging Description

Part Number (Commercial Unit)	F174401
Tip qty/Rack	Rack of 96 tips
Tip qty/Commercial Unit (CU)	10 racks of 96 tips (960 tips)
Tip qty/Outer Package	5 CU of 960 tips (4800 tips)
Codification Lot #	YYMMDD - xxxx (DD: Year, MM: Month, DD: Day, - , xxxx: Mould No.)
Color Code Rack	Blue rack