



AmpliPur Expert Tips

Part Number	F174501
-------------	---------

Specification

Description	AmpliPur Expert Tips, 100-1200 µL, graduated, filter tips
Dimensions	<ul style="list-style-type: none">Tip length: 101.5 mmInner diameter tip point: 0.9 mmOutside tip collar diameter: 8.7 mmRack box dimensions : 122.5 x 90.5 x 119 mm
Volume Range	100-1200 µL
Nominal Volume of Compatible Pipettes	1200 µL
Material (Plastic)	<ul style="list-style-type: none">Tips: Polypropylene (PP)Filter: Polyethylene (PE)Rack: Polypropylene (PP)Heat shrink: Plastic Optical Fibre (POF)
Disposable	<ul style="list-style-type: none">Tips: yesRack: reusable
Tip Color	Translucent
Sterilization	N/A
Quality Control	Material control Product inspection: control of characteristics according Gilson's specifications

Characteristics

Characteristic	The optimized geometry and material (virgin polypropylene) ensures full recovery of the aspirated liquid for consistent results
Temperature	Maximum temperature 140° C
Autoclavability	Not applicable on filter tip
Expiration Date/Shelf Life	3 years
Traceability/Quality	<ul style="list-style-type: none">Complete traceability from production to the userCertificate of conformityQuality statement with protocols and detection thresholds for each batch
Other Information	<ul style="list-style-type: none">Compliant to ISO 8655-2 maximum permissible errorsFree of RNase and DNaseCompliant to the use of hazardous substances directive 2015/863/EUGraduation mark at 50 µL, 100 µL, 200 µL, 500 µL, and 1000 µL

Packaging Description

Part Number (Commercial Unit)	F174501
Tip qty/Rack	Rack of 96 tips
Tip qty/Commercial Unit (CU)	10 racks of 96 tips (960 tips)
Tip qty/Outer Package	5 CU of 960 tips (4800 tips)
Codification Lot #	YYMMDD - xxxx (DD: Year, MM: Month, DD: Day, - , xxxx: Mould No.)
Color Code Rack	Blue rack

gilson.com/contactus



LT800611/A | ©2021 Gilson S.A.S.
All rights reserved. All trademarks
are the property of Gilson, Inc.
and its subsidiaries.