

Standard Analysis Certificate

PRODUCT :ACETYLACETONE RPE
CODE :400300
METHOD :9165

TEST	U.M.	SPECIFICATION
Description	-	Clear liquid
Identification (I.R.)	-	Conform
Water miscibility	-	Conform
Density at 20° C	-	0.971 ÷ 0.981
Refractive index at 20°C	-	1.4510 ÷ 1.4540
Water (K.F.)	%	<=0.1
Residue on evaporation	ppm	<=100
Assay (GLC)	%	>=99.5

DOMINIQUE DUTSCHER SAS

Date

:11/01/2020

QUALITY CONTROL RESPONSIBLE**B. COULANGE (VDR)**