

Standard Analysis Certificate

PRODUCT :PROPIONIC ACID RPE
CODE :409550
METHOD :7662

TEST	U.M.	SPECIFICATION
Description	-	Clear colourless liquid
Identification	-	Positive
Water miscibility	-	Conform
Alcohol miscibility	-	Conform
Density at 20° C	-	0.992 ÷ 0.994
Refractive index at 20°C	-	1.3864 ÷ 1.3884
Boiling point	° C	140 ÷ 142
Water (K.F.)	%	<=0.25
Chloride	ppm	<=5
Heavy metals (Pb)	ppm	<=5
Residue on evaporation	ppm	<=50
Sulphate	ppm	<=10
Fe	ppm	<=2
Assay (acidimetric)	%	>=99.5

Date

:11/01/2020

QUALITY CONTROL RESPONSIBLE

B. COULANGE (VDR)