

# TECHNICAL DATA SHEET

---

<b>Part Number</b>	40967-10X1ML-HON
<b>Description</b>	Trifluoroacetic acid

---

## Technical specifications :

- ✓ trifluoroacetic acid, eluent additive for lc-ms 1ml x10
- ✓ CAS N°76-05-1
- ✓ Molecular weight 114.02 g/mol
- ✓ Molecular formula CF<sub>3</sub>COOH
- ✓ Density 1.48 g/cm<sup>3</sup> (20 °C)
- ✓ Flash point 100 °C
- ✓ Material form liquid
  
- ✓ Storage condition Room Temperature
- ✓ Package type Ampoule
- ✓ Package size 10 x 1mL
  
- ✓ GHS pictograms GHS07 ; GHS05
- ✓ H Statements H412 ; H314 ; H332
- ✓ P Statements P284 ; P280 ; P260

## Sales unit :

10 x 1mL Trifluoroacetic acid

---