

Standard Analysis Certificate

PRODUCT :SULPHURIC ACID 1 N (0.5 mol/L) RPE**CODE** :410570000

TEST	U.M.	SPECIFICATION
Description Assay (potentiometry)	- N	Clear colourless liquid 0.999 - 1.001

his volumetric solution is checked by Tris(hydroxyméthyl)-aminomethan volumetric standard, traceable to the SRM 723e of NIST (National Institute of Standards and Technology)

DOMINIQUE DUTSCHER SAS

Date

:11/01/2020

QUALITY CONTROL RESPONSIBLE**B. COULANGE (VDR)**