

New Era Pump Systems, Inc.

Tel 631-249-1392
Fax 707-248-2089

138 Toledo Street
Farmingdale, NY
11735-6625

www.syringepump.com
info@syringepump.com



CE-DECLARATION OF CONFORMITY

We affirm that the electrical equipment manufactured by us fulfils the requirements of the following EC directives where applicable:

Application of Council Directive(s)	2006/95/EC Low Voltage Directive (LVD) 2004/108/EC EMC Directive
Referenced standards or specifications	Safety: EN61010-1 EMC: EN61236-1:2013 IEC 61000-4-2:2008 / EN61000-4-2:2009 IEC 61000-4-3:2006 / EN61000-4-3:2006 IEC 61000-4-4:2012 / EN61000-4-4:2012 IEC 61000-4-5:2005 / EN61000-4-5:2006 IEC 61000-4-6:2008 / EN61000-4-6:2009 IEC 61000-4-11:2004 / EN61000-4-11:2004 EN55011:2009 / AI:2010 Group 1 Class A ICES-003-Issue 2:2009 Class A FCC Part 15 Class A ANCI C63.4:2014

Model Number	Equipment type	Year of Manufacture
NE-300	Laboratory Infusion Pump	2016

Manufacturer:
NEW ERA PUMP SYSTEMS, INC.
138 Toledo St. Farmingdale, NY 11735 U.S.A


I, the undersigned, hereby declare to my best ability that the equipment specified above conforms to the above Directive(s)

Farmingdale, NY

Place

7/28/16

Date


Signature

Barry Cowen

Full Name

President

Position