

Standard Analysis Certificate

PRODUCT :SULPHURIC ACID 0.5 N (0.25 mol/L) RPE
CODE :410660000

TEST	U.M.	SPECIFICATION
Description Assay (potentiometry)	- N	Clear colourless liquid 0.4995 - 0.5005

This volumetric solution is checked by Tris(hydroxyméthyl)-aminomethan volumetric standard, traceable to the SRM 723e of NIST (National Institute of Standards and Technology)

DOMINIQUE DUTSCHER SAS

Date :11/01/2020

QUALITY CONTROL RESPONSIBLE
B. COULANGE (VDR)