

# Standard Analysis Certificate

**PRODUCT** :SULPHURIC ACID 0.1 N (0.05 mol/L) RPE**CODE** :410710000

TEST	U.M.	SPECIFICATION
Description Assay (potentiometry)	- N	Clear colourless liquid 0.0999 - 0.1001

his volumetric solution is checked by Tris(hydroxyméthyl)-aminomethan volumetric standard, traceable to the SRM 723e of NIST (National Institute of Standards and Technology)

DOMINIQUE DUTSCHER SAS

**Date**

:11/01/2020

**QUALITY CONTROL RESPONSIBLE****B. COULANGE (VDR)**