

5 μ L, Neuros Model 75 RN, point style 3, SYR

Specification Sheet (Part# 65460-02)



Product Specifications

Volume	5 μ L
Syringe Type	Microliter
Syringe Series	700
Termination	Neuros
Needle	Small RN
Gauge	33 gauge
Point Style	3
Needle Length	0 - 20 mm adjustable
Needle Inner Diameter	0.004 inches (0.108 mm)
Needle Outer Diameter	0.008 inches (0.210 mm)
Std. Needle Dead Volume	0.79 μ L
Fluid Path	SST, borosilicate glass
Autoclavable	Yes
Gas Sterilizable	Yes
Pressure	2000 psig (137.9 bar)
Maximum Temperature	115°C / 239°F

Product Specifications

Minimum Temperature	10°C / 50°F
Barrel Inner Diameter	0.014 inches (0.343 mm)
Barrel Outer Diameter	0.260 inches (6.60 mm)
Sales Quantity	1 syringe per pack

HAMILTON 

Web: www.hamiltoncompany.com

USA: +1-800-648-5950

Europe: +41 58 610 10 10

Hamilton Americas & Pacific Rim

Hamilton Company Inc.
4970 Energy Way
Reno, Nevada 89502 U.S.A.
Tel: +1-800-648-5950
Fax: +1-775-856-7259
sales@hamiltoncompany.com

Hamilton Europe, Asia, & Africa

Hamilton Bonaduz A.G.
Via Crusch 8
CH-7402 Bonaduz, Switzerland
Tel: +41 58 610 10 10
Fax: +41-58-610-00-10
contact@hamilton.ch

©2017 Hamilton Company. All Rights Reserved.
All trademarks are owned and/or registered by Hamilton Company in the U.S. and/or other countries.

Revision: A / 2010-05-11 Printed in the U.S.A.

To find a representative in your area, please visit hamiltoncompany.com/contacts