

# Standard Analysis Certificate

**PRODUCT** :ACETONE RS For HPLC  
**CODE** :412500  
**METHOD** :8297

TEST	U.M.	SPECIFICATION
Description	-	Clear colourless liquid
Identification	-	Positive
Density at 20° C	-	0.790 ÷ 0.792
Refractive index at 20°C	-	1.3581 ÷ 1.3601
Boiling point	° C	55.8 ÷ 56.3
Water (K.F.)	ppm	<= 500
Residue on evaporation	ppm	<=5
Acidity	meq/g	<=0.0005
Alcalinity	meq/g	<=0.0002
Assay (GLC)	%	>=99.9
U.V. Transmittance	-	-
At 335 nm	%	>= 60
at 340 nm	%	>= 85
At 345 nm	%	>= 90
at 350 nm	%	>= 98
at 360 nm	%	>= 99

Date

:11/01/2020

QUALITY CONTROL RESPONSIBLE

B. COULANGE (VDR)