

# Standard Analysis Certificate

**PRODUCT** :ETHANOL ANHYDRE RS For HPLC  
**CODE** :412520  
**METHOD** :8298

TEST	U.M.	SPECIFICATION
Description	-	Clear colourless liquid
Identification	-	Positive
Density at 20° C	-	0.7893 ÷ 0.7899
Refractive index at 20°C	-	1.3602 ÷ 1.3622
Boiling point	° C	78.3 ÷ 78.8
Acidity or alkalinity	meq/g	<=0.0002
Water (K.F.)	%	<=0.05
Residue on evaporation	ppm	<=5
Assay (GLC)	%	>=99.9
U.V. Transmittance	-	-
at 210 nm	%	>=25
at 220 nm	%	>=50
at 230 nm	%	>=75
At 240 nm	%	>= 80
at 250 nm	%	>=90
at 270 nm	%	>=94
at 290 nm	%	>=97

Filtered at 0.2 micron

Date

:11/01/2020

**QUALITY CONTROL RESPONSIBLE**

**B. COULANGE (VDR)**