

Standard Analysis Certificate

PRODUCT :ETHYL ACETATE RS For HPLC
CODE :412610000
METHOD :7250

TEST	U.M.	SPECIFICATION
Description	-	Clear colourless liquid
Identification	-	Positive
Density at 20° C	-	0.898 ÷ 0.902
Refractive index at 20°C	-	1.3699 ÷ 1.3739
Boiling point	° C	76.9 ÷ 77.4
Water (K.F.)	ppm	<=300
Residue on evaporation	ppm	<=10
Acidity or alkalinity	meq/g	<=0.0015
Assay (GLC)	%	>=99.9
U.V. Transmittance	-	-
at 260 nm	%	>=76
at 270 nm	%	>=94
at 300 nm	%	>= 97
Methyl alcohol	ppm	<= 100
Ethyl alcohol	ppm	<= 400

Filtered at 0.2 micron

DOMINIQUE DUTSCHER SAS

Date

:11/01/2020

QUALITY CONTROL RESPONSIBLE

B. COULANGE (VDR)