

Standard Analysis Certificate

PRODUCT :DICHLOROMETHANE RS For HPLC
CODE :412620000
METHOD :7248

TEST	U.M.	SPECIFICATION
Description	-	Clear colourless liquid
Identification	-	Positive
Density at 20° C	-	1.322 ÷ 1.328
Refractive index at 20°C	-	1.4214 ÷ 1.4274
Boiling point	° C	39.6 ÷ 40.1
Acidity	meq/g	<= 0.0001
Water (K.F.)	ppm	<=100
Residue on evaporation	ppm	<=2
Assay (GLC)	%	>=99.9
U.V. Transmittance	-	-
At 235 nm	%	>= 40
at 240 nm	%	>=70
at 250 nm	%	>=96
at 260 nm	%	>=99
Amylene	ppm	<= 60
Alkalinity	meq/g	<= 0.0002

Filtered at 0.2 micron

Date

:11/01/2020

QUALITY CONTROL RESPONSIBLE

B. COULANGE (VDR)