

Standard Analysis Certificate

PRODUCT :TOLUENE RS FOR HPLC
CODE :412640000
METHOD :8300

TEST	U.M.	SPECIFICATION
Description	-	Clear colourless liquid
Identification	-	Positive
Density at 20° C	-	0.865 ÷ 0.869
Refractive index at 20°C	-	1.4931 ÷ 1.4991
Boiling point	° C	110.1 ÷ 111.1
Acidity or alkalinity	meq/g	<=0.0015
Water (K.F.)	ppm	<=100
Residue on evaporation	ppm	<=2
Assay (GLC)	%	>=99.8
U.V. Transmittance	-	-
At 285 nm	%	>= 10
at 290 nm	%	>=55
at 300 nm	%	>=80
at 310 nm	%	>=90
at 320 nm	%	>=94
at 330 nm	%	>=96
at 350 nm	%	>=99

Filtered at 0.2 micron

DOMINIQUE DUTSCHER SAS

Date

:11/01/2020

QUALITY CONTROL RESPONSIBLE

B. COULANGE (VDR)