

# Standard Analysis Certificate

**PRODUCT** :DICHLOROMEHANE STAB ETHANOL RS HPLC  
**CODE** :412660  
**METHOD** :7278

TEST	U.M.	SPECIFICATION
Description	-	Clear colourless liquid
Identification	-	Positive
Boiling point	°C	39.55 ÷ 40.05
Density at 20°C	-	1.322 ÷ 1.328
Refractive index at 20°C	-	1.4214 ÷ 1.4274
Acidity (Formic acid)	ppm	<= 5
Residue on evaporation	ppm	<= 5
Ethanol	%	0.1 ÷ 0.4
Water (K.F.)	ppm	<= 100
Transmittance U.V.	-	
at 235 nm	%	>= 40
at 240 nm	%	>= 75
at 245 nm	%	>= 90
at 255 nm	%	>= 96
at 260 nm	%	>= 99
Assay (GLC)	%	>= 99.9
Filtered at 0.2 µm	-	

Date

:11/01/2020

QUALITY CONTROL RESPONSIBLE

B. COULANGE (VDR)