

# Standard Analysis Certificate

**PRODUCT** :DIETHYL ETHER RS HPLC  
**CODE** :412670  
**METHOD** :412670

TEST	U.M.	SPECIFICATION
Clear, colourless liq. appearance	-	Conform
Identification	-	Conform
Colour	Apha	<= 10
Refractive index at 20°C	-	1.35 - 1.354
Water content (K.F.)	mg/Kg	<= 100
Free acid (as CH <sub>3</sub> COOH)	mg/Kg	<= 30
Peroxides (as H <sub>2</sub> O <sub>2</sub> )	mg/Kg	<= 60
Non volatile residue	mg/Kg	<= 5
UV transmittance at 220 nm	%	>= 10
UV transmittance at 230 nm	%	>= 45
UV transmittance at 250 nm	%	>= 75
UV transmittance at 280 nm	%	>= 92
UV transmittance at 300 nm	%	>= 98
Assay (GC)	%	>= 99.7

Not stabilized.

DOMINIQUE DUTSCHER SAS

Date :11/01/2020

**QUALITY CONTROL RESPONSIBLE**  
**B. COULANGE (VDR)**