

Standard Analysis Certificate

PRODUCT :2-METHYLTETRAHYDROFURAN HPLC
CODE :412680

TEST	U.M.	SPECIFICATION
Clear, colourless liq. appearance	-	Conform
Identification	-	Conform
Colour	Apha	<= 10
Refractive index at 20°C	-	1.404 - 1.408
Water content (K.F.)	mg/Kg	<= 200
Non volatile residue	mg/Kg	<= 5
Peroxides (as H ₂ O ₂)	mg/Kg	<= 300
Assay (GC)	%	>= 99.5
UV transmittance at 240 nm	%	>= 30
UV transmittance at 250 nm	%	>= 50
UV transmittance at 260 nm	%	>= 70
UV transmittance at 280 nm	%	>= 90
UV transmittance from 310 nm	%	>= 98

Not stabilized.

Date

:11/01/2020

QUALITY CONTROL RESPONSIBLE
B. COULANGE (VDR)