

Standard Analysis Certificate

PRODUCT :BENZYL ALCOHOL RPE
CODE :414020
METHOD :8047

TEST	U.M.	SPECIFICATION
Description	-	Clear liquid
Colour (APHA)	-	<=10
Identification	-	Positive
Density at 20° C	-	~ 1.04
Refractive index at 20° C	-	~ 1.54
Boiling point	°C	204.5 ÷ 205.5
Water (K.F.)	%	<=0.1
Total chlorine	ppm	<=100
Benzaldehyde	%	<=0.1
Assay (GLC)	%	>=99.8

Temperature storage : between 2°C and 8°C

Date

:11/01/2020

QUALITY CONTROL RESPONSIBLE
B. COULANGE (VDR)