

Standard Analysis Certificate

PRODUCT :BUTANOL-2 RPE - Reag. Ph.Eur.
CODE :414260
METHOD :6490

| TEST | U.M. | SPECIFICATION |
|--------------------------|------|-----------------|
| Description | - | Clear liquid |
| Colour (APHA) | - | <=10 |
| Identification | - | Positive |
| Water solubility | - | Conform |
| Density at 20° C | - | 0.801 ÷ 0.811 |
| Refractive index at 20°C | - | 1.3944 ÷ 1.3984 |
| Boiling point | °C | 99 - 100 |
| Water (K.F.) | % | <=0.2 |
| Residue on evaporation | ppm | <=20 |
| Acidity (butirric acid) | ppm | <=20 |
| Alcalinity (NaOH) | ppm | <=10 |
| Indole base | ppm | <=0.1 |
| Assay (GLC) | % | >=99.0 |

Date

:11/01/2020

QUALITY CONTROL RESPONSIBLE

B. COULANGE (VDR)