

Standard Analysis Certificate

PRODUCT :PROPANOL-1 RPE - Reag. Ph.Eur.
CODE :415100
METHOD :6438

TEST	U.M.	SPECIFICATION
Description	-	Clear liquid
Identification (I.R.)	-	Conform
Colour	APHA	<= 10
Density at 20° C	-	0.803 ÷ 0.805
Refractive index at 20°C	-	1.3840 ÷ 1.3860
Boiling point	°C	96 ÷ 98
Distillation range	-	96 ÷ 99°C
Acidity (acetic acid)	%	<= 0.03
Water (K.F.)	ppm	<= 1000
Residue on evaporation	ppm	<= 5
Assay (GLC)	%	>= 99.5

DOMINIQUE DUTSCHER SAS

Date

:11/01/2020

QUALITY CONTROL RESPONSIBLE

B. COULANGE (VDR)