

Technical Datasheet

Article No. 418

Water for HPLC

For laboratory use.

Parameter	Value
Appearance	colourless liquid
Density	1.00 g/ml
Mol.-Weight	18.02 g/mol
Melting point	0 °C
Boiling point	100 °C
HPLC gradient test	conform
HPLC gradient (peak) at 210 nm	max. 5 mAU
HPLC gradient (peak) at 254 nm	max. 0.5 mAU
Conductivity	max. 0.1 µS/cm
Total organic carbon (TOC)	max. 3 ppb
Filtered through 0.1 µm	