

Standard Analysis Certificate

PRODUCT :SILVER NITRATE 0.1 N (0.1 mol/L) RPE**CODE** :424060000

TEST	U.M.	SPECIFICATION
Description Assay (potentiometry)	- N	Clear colourless liquid 0.0999 - 0.1001

This volumetric solution is checked by sodium chloride volumetric standard, traceable to the SRM 399b of NIST (National Institute of Standards and Technology)

Date

:11/01/2020

QUALITY CONTROL RESPONSIBLE**B. COULANGE (VDR)**