



# CERTOCLAV

**Article Name**  
CertoClav Classic

**Article number**  
85.101.74

**Short description**  
Autoclave for sterilisation at 125°C or 140°C. Venting by stream process. Manual process timing with integrated timer. Process control by thermostat and precision valve. 18 liter volume.

### Description

CertoClav EL - the trusted, economic autoclaving solution since 50 years. Process control via precision valve. Now with stainless steel chamber.

The simple operation allows for effortless work with less handling. Through the CertoLock locking device the laborious turning of several screws is avoided.

4-fold safety system; 50 years experience in autoclave construction.

Selective 125/140° C. 100°C programme for the dissolving of culture media.

Manual time control via integrated timer.

Compact instruments also simply packaged, liquids (culture mediums), bottles and equipment, refuse sterilisation.

### Configuration

Precision thermometer	Steam release cock
Precision pressure gauge	Thermostat and regulating valve
4-fold overpressure security system	Heating integrated

### Technical data

Material	Stainless Steel 304L	Certificates	CE, TUEV
----------	----------------------	--------------	----------

### Dimensions

Volume	18 litre	Height	60 cm
Chamber diameter	26cm	Width	44 cm
Chamber height	35cm	Depth	40 cm
Useable height chamber	30 cm	Weight	10 kg

### Working data

Max. working temp	140°C	Current	230V, 50/60 Hz
Max. working pressure	2,7 Bar (39 psi)	Power	1300 Watt
Rated pressure	4,1 Bar (59 psi)		

### Included

Autoclave
Base plate (internal)
Mains cable

### Recommended accessories

Wire basket diam. 25	13.001.49
----------------------	-----------