

# CertoClav MultiControl 2

## THE FREELY PROGRAMMABLE AUTOCLAVE

PRINTER PORT  
MULTI-LANGUAGE  
PC PORT



Art.-No. 8510175 (18 litres)

### FEATURES

- Set programs for all sterilisation tasks
- Logging of processes with PC Software (basic version included) or via printer.
- Simple program programmable on the controller (1 time / temperature step)
- Complex programs (multiple time / temperature steps) can be entered with PC Software (Professional version required).
- Automatic venting process for reliable sterilisation.
- Programmable from room temperature to 140° C
- Suitable for materials testing (ICP-TM610, HAST, GLB-AL 631)

### INCLUDED

- Controller with with printer port and USB port
- Large LCD screen with temperature display

### ACCESSORY

#### Wire basket

Art.-No. 8583010  
Diam. 23, stainless steel  
For easy filling with  
sterilisation goods.



#### Printer

Art.-No. 1300187  
For easy print out of  
your documentation.



#### PC Software Professional Version

Art.-No. 8200210 (Art.-No. 8200211  
when purchased in bundle with autoclave)  
For programming of multi step  
programs.