

Wire basket with handle

Item number 1300149

Description

The sterilization wire basket made of stainless steel for our vertical, stainless steel autoclaves is suitable for solids, waste in bags, instruments, bottles, liquids, material testing and the like. The basket is stackable.

Technical data

Material: Stainless steel 18/10
Diameter: ø 25 cm
Height: 22,5 cm
Weight: 0,6 kg

Accessory for

CertoClav Classic 125/140°C
CertoClav Classic 115/121°C
CertoClav MultiControl 2
CertoClav MultiControl 2 Touch
CertoClav Connect
CertoClav Cibus

Additional product information

Packaging unit: 1 Piece
Cleaning and Reprocessing: Manual cleaning
Warranty: 2 Years



The illustrations, information and technical data correspond to the information available at the time of printing. Technical changes are reserved. Further information is available at www.certoclav.com or personally from our team.