

Standard Analysis Certificate

PRODUCT :CHLOROFORM Stab. Amylene RPE
CODE :438580
METHOD :438580

TEST	U.M.	SPECIFICATION
Clear, colourless liq. appearance	-	Conform
Identification	-	Conform
Colour	Apha	<= 10
Refractive index at 20°C	-	1.444 - 1.448
Water content (K.F.)	mg/Kg	<= 100
Free acid (as HCl)	mg/Kg	<= 5
Non volatile residue	mg/Kg	<= 10
Stabilizer (Amylene)	mg/Kg	5 - 50
Carbon tetrachloride	mg/Kg	<= 80
Dichloromethane	mg/Kg	<= 50
1,2-dichloroethane	mg/Kg	<= 10
Assay (GC)	%	>= 99.95

DOMINIQUE DUTSCHER SAS

Date

:11/01/2020

QUALITY CONTROL RESPONSIBLE

B. COULANGE (VDR)