

Standard Analysis Certificate

PRODUCT :CHLOROFORME stab Ethanol RPE ACS ISO ReagPhEur ReagUSP
CODE :438610
METHOD :7269

TEST	U.M.	SPECIFICATION
Description	-	Clear liquid
Colour (APHA)	-	<=10
Identification (I.R.)	-	Conform
Phosgene	-	Conform
Ready carbonizable substances	-	Conform
Acetone and aldehydes	ACS	Conform
Acidità e cloruri	ACS	Conform
Suit.for Dithizone col.	ACS	Conform
Density at 20° C	-	1.477 ÷ 1.481
Refractive index at 20°C	-	1.4461 ÷ 1.4491
Boiling point	° C	60.5 ÷ 61.5
Water (K.F.)	ppm	<=100
Residue on evaporation	ppm	<=5
Acidity (HCl)	ppm	<=5
Ethyl alcohol	%	0.6 ÷ 1
Chloride	ppm	<=0.4
Free chlorine	ppm	<=0.1
Carbonyl Compounds (CO)	ppm	<=5
Cu	ppm	<=0.01
Fe	ppm	<=0.1
Pb	ppm	<=0.01
Zn	ppm	<=0.05
Assay (GLC)	%	>=99.9
Stab.÷0.75%ethanol anh.	-	

Date

:11/01/2020

QUALITY CONTROL RESPONSIBLE

B. COULANGE (VDR)