

Standard Analysis Certificate

PRODUCT :CHLOROFORM RS Plus Pesticides
CODE :438680

TEST	U.M.	SPECIFICATION
Refractive index at 20°C	-	1.444 - 1.448
Water content (K.F.)	mg/Kg	<= 100
Non volatile residue	mg/Kg	<= 5
Colour	Hazen	<= 10
Free acid (as HCl)	mg/Kg	<= 5
Assay (GC)	%	>= 99.9
Stabilizer (Amylene)	mg/Kg	5 - 50
Retention time trichlorobenzene to mirex	-	-
GC-ECD.Individual peak (Lindane)	ng/l	<= 3
GC-NPD.Individual peak (Ethylparathion)	ng/l	<= 3
Retention time Atrazin to Coumaphos	-	-

For chlorinated and nitrogenous compounds analysis; extraction of pesticides and weed-killers.

Date

:11/01/2020

QUALITY CONTROL RESPONSIBLE
B. COULANGE (VDR)