

Product Specification

Specification	Product Code	128040M
	Description	30ml Type I Neutral Glass Moulded Universal with Fitted Aluminium Cap and Rubber Liner
	Container Material	Type I Neutral Glass
	Colour	Clear
	Capacity	30ml
	Cap Material	Aluminium
	Cap Colour	Gold
	Cap Liner	Bromobutyl Rubber Wad
	Sterile	No
	Autoclavable	Yes
	Universal Dimensions	28.6mm dia x 83mm high (without cap)
	Cap Dimensions	28R3
	Case Quantity	144
	CE Marked	No



(not actual size)



SciLabware is a company of DWK Life Sciences

Pyrex® is a registered trademark; SciLabware is a registered user.