

# PRODUCT SPECIFICATION/CERTIFICATE OF CONFORMITY

<b>Specification</b>	<b>Product Code</b>	327152PP
	<b>Description</b>	30ml Universal Container with Printed label PP Screw Cap
	<b>Material Base</b>	Polypropylene
	<b>Material Cap</b>	Polypropylene
	<b>Colour</b>	Clear
	<b>Sterility</b>	Sterilised by ETO
	<b>Shelf Life</b>	5 Years from Manufacture

<b>Dimension</b>	<b>Height (mm)</b>	93.5
	<b>Diameter (mm)</b>	30.5
	<b>Packaging Format</b>	400/Case 50/Bag

<b>Compliance</b>	<b>Food safe (EU/FDA)</b>	EU 10/2011 / FDA CFR 21 177
	<b>BSE / TSE free</b>	YES- No tallow from animal sources are used in the manufacture of the finished product.
	<b>REACH Compliant</b>	Compliant according to 1907/2006
	<b>LATEX Free</b>	YES- No Materials containing latex natural rubber are used in the manufacturing, handling, and packaging process.
	<b>CE Marking</b>	Meets the essential requirements of European directive 98/79 EC
	<b>Leak Testing</b>	In Accordance with to BS EN 14254 Annexe D



**DWK Life Sciences Limited**

Unit 4, Riverside 2, Campbell Road Stoke-on-Trent, ST4 4RJ, UK

Registered in England and Wales No.7723574

Tel: +44 (0)1782 444406

Fax: +44 (0)1782 940436

[www.dwk Ltd.com](http://www.dwk Ltd.com)  
[enquiries.uk@dwk.com](mailto:enquiries.uk@dwk.com)

DWK Life Sciences Limited hereby confirms that the product specified above is inspected according to all relevant specifications and procedures laid out in our Quality System, in accordance with our registration to BS EN ISO9001 and BS EN ISO13485.

Signed on behalf of DWK Life Sciences Limited by:

Signed: M.Ashley  
Quality Technician

Date: 18 Sep 2020