



### Reagecon pH 1.00 Buffer Solution at 20°C

Item No. 10105

Reagecon's pH Buffer Solution at 20°C is tested by an ISO17025 Accredited Test Method (INAB Ref: 264T) and NIST traceable, has guaranteed stability throughout its entire shelf life, even after opening the bottle.



### Additional Details

Pack Size 500ml  
Weight 0.62  
Weight Unit Kg  
Tariff Code 3822 190090  
Country of Origin Ireland

### Safety Information

UN Number	1789	Packing Group	III
Class	8	PG Tunnel	E

## DESCRIPTION

Tested by Reagecon's ISO 17025 Accredited Test Method (INAB Ref: 264T)

NIST traceable

# PRODUCT SPECIFICATIONS

pH Buffer Standard, pH Buffer Material

Value	<b>pH 1.00 ± 0.02</b>
Additional Information	<b>Accuracy ± 0.02 pH units</b>
Tested	<b>at 20°C</b>

Use	<b>Laboratory Chemical</b>
Packaging	<b>HDPE Bottle GL45</b>
Physical Form	<b>Liquid</b>

DOMINIQUE DUTSCHER SAS