



Reagecon pH 4.00 Buffer Solution at 20°C

Item No. 10405

Reagecon's pH Buffer Solution at 20°C is tested by an ISO17025 Accredited Test Method (INAB Ref: 264T) and NIST traceable, has guaranteed stability throughout its entire shelf life, even after opening the bottle.



Additional Details

Pack Size	500ml
Weight	0.62
Weight Unit	Kg
Tariff Code	3822 190090
Country of Origin	Ireland

DESCRIPTION

Tested by Reagecon's ISO 17025 Accredited Test Method (INAB Ref: 264T)

NIST traceable

PRODUCT SPECIFICATIONS

pH Buffer Standard, pH Buffer Material

Value	pH 4.00 ± 0.01
Additional Information	Accuracy ± 0.01 pH units
Tested	at 20°C

Use	Laboratory Chemical
Packaging	HDPE Bottle GL45
Physical Form	Liquid

DOMINIQUE DUTSCHER SAS