



Delivering the Correct Result



## Reagecon pH 4.00 Colour Coded Buffer Solution at 20°C

Item No. 1040C

Use Colour Coded pH buffer for ease of use. Reagecon's pH Colour Coded Buffer Solution tested at 20°C is tested by an ISO17025 Accredited Test Method (INAB Ref: 264T) and NIST traceable, has guaranteed stability throughout its entire shelf life, even after opening the bottle.

### Additional Details

Pack Size 1L  
Weight 1.16  
Weight Unit Kg  
Colour Red

Country of Origin Ireland

## DESCRIPTION

Tested by Reagecon's ISO 17025 Accredited Test Method (INAB Ref: 264T)

NIST traceable

DOMINIQUE DUTSCHER SAS

# PRODUCT SPECIFICATIONS

pH Buffer Standard, pH Buffer Material

Value	<b>pH 4.00 ± 0.01</b>
Additional Information	<b>Accuracy ± 0.01 pH units</b>
Tested	<b>at 20°C</b>

Use	<b>Laboratory Chemical</b>
Packaging	<b>HDPE Bottle GL45</b>
Physical Form	<b>Liquid</b>

DOMINIQUE DUTSCHER SAS