



Reagecon pH 7.00 Colour Coded Buffer Solution at 25°C

Item No. 107025C

Use Colour Coded pH buffer for ease of use. Reagecon's colour coded pH Buffer Solution at 25°C is tested and certified by an ISO17025 Accredited Test Method (INAB Ref: 264T) and NIST traceable, has guaranteed stability throughout its entire shelf life, even after opening the bottle.

Additional Details

Pack Size	1L
Weight	1.16
Weight Unit	Kg
Colour	Yellow
Tariff Code	3822 0000
Country of Origin	Ireland

DESCRIPTION

- Certified by Reagecon's ISO 17025 Accredited Test Method (INAB Ref: 264T)
- NIST traceable

PRODUCT SPECIFICATIONS

- pH Buffer Standard, pH Buffer Material

Value	pH 7.00 ± 0.01
Additional Information	Accuracy ± 0.01 pH units
Certified	at 25°C

Use	Laboratory Chemical
Packaging	HDPE Bottle GL45
Physical Form	Liquid