



## Reagecon pH 10.00 Colour Coded Buffer Solution at 25°C

Item No. 110025C

Use Colour Coded pH buffer for ease of use. Reagecon's colour coded pH Buffer Solution at 25°C is tested and certified by an ISO17025 Accredited Test Method (INAB Ref: 264T) and NIST traceable, has guaranteed stability throughout its entire shelf life, even after opening the bottle.

### Additional Details

Pack Size	1L
Weight	1.16
Weight Unit	Kg
Colour	Blue
Tariff Code	3822 0000
Country of Origin	Ireland

## DESCRIPTION

- Certified by Reagecon's ISO 17025 Accredited Test Method (INAB Ref: 264T)
- NIST traceable

## PRODUCT SPECIFICATIONS

- pH Buffer Standard, pH Buffer Material

Value	<b>pH 10.00 ± 0.01</b>
Additional Information	<b>Accuracy ± 0.01 pH units</b>
Certified	<b>at 25°C</b>

Use	<b>Laboratory Chemical</b>
Packaging	<b>HDPE Bottle GL45</b>
Physical Form	<b>Liquid</b>