

Reagecon

A CALIBRE SCIENTIFIC COMPANY



Reagecon Osmolality Standard 500 mOsm/Kg

Item No. RE-OSM-500



Reagecon's Osmolality Standard 500 mOsm/Kg is used to calibrate, control or validate all types of Osmometers. The Osmolality results are reported in milliOsmoles per kilogram water (mOsm/kg H₂O shortened to mOsm/kg for routine use). This product is manufactured in compliance with European, United States and Chinese Pharmacopoeia guidelines where appropriate.

Additional Details

Pack Size	12x5ml
Weight	0.72
Weight Unit	Kg
Tariff Code	3822 190090
Country of Origin	Ireland

DESCRIPTION

Tested by Reagecon's ISO 17025 Accredited Test Method (INAB Ref: 264T)

This Osmolality Standard is prepared gravimetrically on a weight/weight basis. Both solute and solvent are weighed on a balance calibrated by Reagecon engineers using OIML traceable weights. Reagecon holds ISO 17025 accreditation for calibration of laboratory balances (INAB Ref:265C).

PRODUCT SPECIFICATIONS

Osmolality Standard, Osmolality Material

Value	500 mOsm/Kg
Packaging	Amber Glass Ampoule
Physical Form	Liquid

Use	Osmometer Calibration
Tested	Yes

DOMINIQUE DUTSCHER SAS