

# Reagecon

A CALIBRE SCIENTIFIC COMPANY



## Reagecon Osmolality Standard 600 mOsm/Kg

Item No. RE-OSM-600



Reagecon's Osmolality Standard 600 mOsm/Kg is used to calibrate, control or validate all types of Osmometers. The Osmolality results are reported in milliOsmoles per kilogram water (mOsm/kg H<sub>2</sub>O shortened to mOsm/kg for routine use). This product is manufactured in compliance with European, United States and Chinese Pharmacopoeia guidelines where appropriate.

### Additional Details

Pack Size	12x5ml
Weight	0.72
Weight Unit	Kg
Tariff Code	3822 190090
Country of Origin	Ireland

## DESCRIPTION

Tested by Reagecon's ISO 17025 Accredited Test Method (INAB Ref: 264T)

This Osmolality Standard is prepared gravimetrically on a weight/weight basis. Both solute and solvent are weighed on a balance calibrated by Reagecon engineers using OIML traceable weights. Reagecon holds ISO 17025 accreditation for calibration of laboratory balances (INAB Ref:265C).

# PRODUCT SPECIFICATIONS

Osmolality Standard, Osmolality Material

Value	<b>600 mOsm/Kg</b>
Packaging	<b>Amber Glass Ampoule</b>
Physical Form	<b>Liquid</b>

Use	<b>Osmometer Calibration</b>
Tested	<b>Yes</b>

DOMINIQUE DUTSCHER SAS