

This website uses cookies to enhance your experience. [I agree](#) [More information](#)

# Reagecon

## Delivering the Correct Result



### Reagecon Cryoscope Heat Transfer Fluid

Item No. HTF250

Reagecon's Cryoscope Heat Transfer Fluid is used as the temperature control liquid in Cryoscope cooling systems. The product is manufactured for use on all Cryoscopes compliant to International Reference Standard ISO5764/IDF108 for the determination of freezing point in milk.

#### Additional Details

Pack Size	250ml
Weight	0.30
Weight Unit	Kg
Tariff Code	3822 0000
Country of Origin	Ireland

## DESCRIPTION

The Heat Transfer Fluid is formulated to ensure the transfer of thermal energy and completion of the freezing point process accurately. It is presented in 250ml bottles designed to be used in conjunction with specific Cryoscope types.

# PRODUCT SPECIFICATIONS

Cryoscope Standard, Cryoscope Material

Use	<b>Cryoscopy</b>
Tested	<b>Yes</b>

Packaging	<b>HDPE Bottle</b>
Physical Form	<b>Liquid</b>

DOMINIQUE DUTSCHER SAS