

Reagecon

A CALIBRE SCIENTIFIC COMPANY



Reagecon 100 ppm Chemical Oxygen Demand Standard

Item No. COD100



This highly accurate 100 ppm Chemical Oxygen Demand (COD) Calibration Standard, can be used to calibrate, control, verify or validate all COD Test Methods. The product yields an accurate test result of 100 ppm. Achieving this correct result increases significantly the confidence of the analyst in the correct COD reading of the sample.

Additional Details

Pack Size	500ml
Weight	0.62
Weight Unit	Kg
Tariff Code	3822 190090
Country of Origin	Ireland

DOMINIQUE DUTSCHER SAS

PRODUCT SPECIFICATIONS

Chemical Oxygen Demand Standard, Chemical Oxygen Demand Material

Value	100 ppm COD
Packaging	HDPE Bottle GL45
Physical Form	Liquid

Use	For COD Test Methods
Tested	Yes

DOMINIQUE DUTSCHER SAS