

# Reagecon

A CALIBRE SCIENTIFIC COMPANY



## Reagecon 500 ppm Chemical Oxygen Demand Standard

Item No. COD500



This highly accurate 500 ppm Chemical Oxygen Demand (COD) Calibration Standard, can be used to calibrate, control, verify or validate all COD Test Methods. The product yields an accurate test result of 500 ppm. Achieving this correct result increases significantly the confidence of the analyst in the correct COD reading of the sample.

### Additional Details

Pack Size	500ml
Weight	0.62
Weight Unit	Kg
Tariff Code	3822 190090
Country of Origin	Ireland

DOMINIQUE DUTSCHER SAS

# PRODUCT SPECIFICATIONS

Chemical Oxygen Demand Standard, Chemical Oxygen Demand Material

Value	<b>500 ppm COD</b>
Packaging	<b>HDPE Bottle GL45</b>
Physical Form	<b>Liquid</b>

Use	<b>For COD Test Methods</b>
Tested	<b>Yes</b>

DOMINIQUE DUTSCHER SAS