

Standard Analysis Certificate

PRODUCT :DICHLOROMETHANE Stab.amylene RS For Spectroscopy
CODE :442370
METHOD :8503

TEST	U.M.	SPECIFICATION
Description	-	Clear colourless liquid
Colour (APHA)	-	<=10
Identification (I.R.)	-	Positive
Density at 20° C	-	1.322 ÷ 1.328
Refractive index at 20°C	-	1.4214 ÷ 1.4274
Boiling point	° C	39.6 ÷ 40.1
Water (K.F.)	ppm	<=100
Residue on evaporation	ppm	<=5
Acidity	meq/g	<=0.0005
Assay (GLC)	%	>=99.9
Fluorescence	-	-
at 254 nm	ppb	<=2
at 365 nm	ppb	<=2
U.V. Transmittance	-	-
at 235 nm	%	>=35
at 240 nm	%	>=70
at 245 nm	%	>=90
at 250 nm	%	>=95
at 260 nm	%	>=99
Amylene	ppm	20 ÷ 60
Alkalinity	meq/g	<= 0.0002

DOMINIQUE DUTSCHER SAS

Date :11/01/2020

QUALITY CONTROL RESPONSIBLE
B. COULANGE (VDR)