

TECHNICAL DATASHEET

REFERENCE :

ClearLine® Test Strips

The ClearLine® test strips are designed for semi-quantitative analysis in water. Most tests can be performed under 90 seconds with these strips, providing a fast and accurate reading every time.

They are ready-to-use, pre-calibrated and contain all necessary equipment and reagents. The colour schemes associated with different tests are adjusted and checked using certified standard solutions, traceable to industry standards.



Ordering Information :

Product	Description	Sales Unit (Strips)	Cat No
Ammonium 400	0 - 10 - 25 - 50 - 100 - 200 - 400 mg/L NH ₄ ⁺	50	442414
Ascorbic Acid 2000	0 - 50 - 100 - 200 - 300 - 500 - 700 - 1000 - 2000 mg/L	100	442415
Peracetic Acid 50	0 - 5 - 10 - 20 - 30 - 50 mg/L	100	442409
Peracetic Acid 500	0 - 100 - 150 - 200 - 250 - 300 - 400 - 500 mg/L	100	442410
Peracetic Acid 2000	0 - 50 - 1000 - 1500 - 2000 mg/L	100	442422
Chlorine 5	0 - 0,5 - 1 - 3 - 5 mg/L Cl ₂	100	442405
Chlorine 20	0 - 1 - 3 - 5 - 10 - 20 mg/L Cl ₂	100	442416
Chlorine 300	0 - 25 - 50 - 100 - 200 - 300 mg/L Cl ₂	100	442406
Chlorine 1000	0 - 50 - 100 - 250 - 500 - 1000 mg/L	100	442407
Water Hardness	3 - 4 - 7 - 14 - 21°d	100	442421
Nitrate 500	0 - 10 - 25 - 50 - 250 - 500 mg/L	100	442413
Nitrite 25	0 - 0,5 - 1 - 5 - 10 - 25 mg/L	100	442412
Nitrite 80	0 - 2 - 5 - 10 - 20 - 40 - 80 mg/L	100	442417
Peroxyde 25	0 - 0,5 - 2 - 5 - 10 - 25 mg/L	100	442418
Peroxyde 100	0 - 1 - 3 - 10 - 30 - 100 mg/L	100	442408
Peroxyde 1000	0 - 100 - 200 - 400 - 600 - 800 - 1000 mg/L	100	442419
Phosphate 100	0 - 3 - 10 - 25 - 50 - 100 mg/L	50	442420
QAC 400	0 - 100 - 200 - 300 - 400 mg/L	100	442404
Sulfite 500	10 - 50 - 100 - 250 - 500 mg/L	100	442411