

Standard Analysis Certificate

PRODUCT :1,4-DIOXANE RS HPLC Isocratique- Stabilized
CODE :443230
METHOD :443230

| TEST | U.M. | SPECIFICATION |
|------------------------|------|---------------|
| Description | - | Clear liquid |
| Colour | APHA | <= 20 |
| Identification (I.R.) | - | Positive |
| Water (K.F.) | ppm | <= 250 |
| Residue on evaporation | ppm | <= 5 |
| Acidity (acetic acid) | ppm | <= 50 |
| Alcalinity (NH3) | ppm | <= 3 |
| Assay (CPG) | % | >= 99.5 |
| Absorbance | - | - |
| At 245 nm | AU | <= 0.229 |
| At 270 nm | AU | <= 0.097 |
| Transmittance | - | - |
| At 245 nm | % | >= 59 |
| At 270 nm | % | >= 80 |

Stab.with 1 ppm BHT

DOMINIQUE DUTSCHER SAS

Date :11/01/2020

QUALITY CONTROL RESPONSIBLE
B. COULANGE (VDR)